

## Automated External Defibrillator with Video Display

By Staff Editor

Aug 26, 2010 - 7:43:53 PM



(HealthNewsDigest.com) - The Defibtech Lifeline View™ Automated External Defibrillator (AED) is the newest addition to the Lifeline AED family of products. It is the first and only AED available with a full-color, full-motion interactive video display that shows step-by-step videos for performing CPR, rescue breathing and external defibrillation. It makes it easy for anyone—from first responders to untrained bystanders—to respond confidently and appropriately during an emergency. The device combines a lightweight, user-friendly design and weighs less than three pounds and is tested to U.S. Military “drop and shock” specifications. The interactive video screen shows the user what to do with clear, easy to understand video, pictures and text and a calm voice so the user can see it, read it and hear it. This means that it’s easy to know what to do to help save a life, from preparing a patient and using the AED correctly, to performing CPR and rescue breathing.

For more information visit: <http://defibview.com/>