

# Lifeline ARM

## Automated Chest Compression System for Professionals

With an innovative and elegant design, Defibtech's latest technological innovation provides victims of sudden cardiac arrest the precision needed for high-quality and continuous CPR associated with better survival. Easy to deploy and use, the device delivers compressions, with complete chest recoil, at the depth and rate as recommended by AHA/ERC cardiopulmonary resuscitation (CPR) guidelines. Using a proprietary algorithm, the removable compression module delivers smooth and consistent compressions. A rigid backboard with a stiff frame of single-piece design facilitates operation without unwanted flex.

***Be it on the ground, in an ambulance cot, a moving vehicle, or intra-hospital transport, the Lifeline ARM is your solution for uninterrupted CPR.***



- Removable Compression Module allows for quick deployment, enhanced usability
- Built for the demanding environments of the professional first responder
- Lightweight rigid frame
- Compensates for variability in chest resistance
- Long-life rechargeable battery pack



Defibtech, LLC | 741 Boston Post Road, Guilford, CT 06437 USA  
www.defibtech.com | 1-866-DEFIB-4U (1-866-333-4248) or 1-203-453-4507

**USA** Rx ONLY

RAC-A1703EN-BB  
**ELECTRONIC DISTRIBUTION**

THIS DOCUMENT IS FOR USA DISTRIBUTION ONLY

The ARM is intended for use by qualified medical personnel, certified to administer CPR, as an adjunct to manual CPR when effective manual CPR is not possible (e.g., during patient transport, or extended CPR when fatigue may prohibit the delivery of effective/consistent compressions to the victim, or when insufficient personnel are available to provide effective CPR).