

Lifeline ARM

Automated Chest Compression System for Professionals

With an innovative and elegant design, Defibtech's latest technological innovation provides victims of sudden cardiac arrest the precision needed for high-quality and continuous CPR associated with better survival. Easy to deploy and use, the device delivers compressions, with complete chest recoil, at the depth and rate as recommended by AHA/ERC cardiopulmonary resuscitation (CPR) guidelines. Using a proprietary algorithm, the removable compression module delivers smooth and consistent compressions. A rigid backboard with a stiff frame of single-piece design facilitates operation without unwanted flex.



Be it on the ground, in an ambulance cot, a moving vehicle, or intra-hospital transport, the Lifeline ARM is your solution for uninterrupted CPR.



- Removable Compression Module allows for quick deployment, enhanced usability
- Built for the demanding environments of the professional first responder
- Lightweight rigid frame
- Compensates for variability in chest resistance
- Long-life rechargeable battery pack



Defibtech, LLC | 741 Boston Post Road, Guilford, CT 06437 USA
www.defibtech.com | 1-866-DEFIB-4U (1-866-333-4248) or 1-203-453-4507

RAC-E1703EN-BB

ELECTRONIC DISTRIBUTION