

IP Ratings of Defibtech Lifeline AEDs

- Defibtech AEDs are highly rated for protection against dust and water ingress
- Defibtech 100-Series Lifeline AEDs have an IP rating of **IP54**
- Defibtech 2000-Series Lifeline AEDs all have an IP rating of **IP55**

NEMA Standard IEC 60529 IP Ratings (IP Ratings standard for dust exposure and water intrusion).

| IP First number – Protection against solid objects | IP Second number – Protection against liquids |
|--|--|
| 0 Non-protected | 0 Non-protected |
| 1 Protected against solid foreign objects of 50 mm diameter and greater | 1 Protected against vertically falling water drops |
| 2 Protected against solid foreign objects of 12.5 mm diameter and greater | 2 Protection against vertically falling water drops when enclosure tilted up to 15° |
| 3 Protected against solid foreign objects of 2.5 mm diameter and greater | 3 Protected against spraying water |
| 4 Protected against solid foreign objects of 1.0 mm diameter and greater | 4 Protected against splashing water |
| 5 Dust protected | 5 Protected against water jets |
| 6 Dust tight | 6 Protected against powerful water jets |
| | 7 Protected against the effects of temporary immersion in water |
| | 8 Protects against the effects of continuous immersion in water |

Ingress Protection (IP) ratings are developed by the European Committee for Electro Technical Standardization (CENELEC) (NEMA IEC 60529 Degrees of Protection Provided by Enclosures - IP Code), specifying the environmental protection the enclosure provides.



Life-Saving Technology Within Reach

Defibtech, LLC • 741 Boston Post Road, Guilford, CT 06437 USA
www.defibtech.com • 1-866-DEFIB-4U • 1-203-453-4507

DAC-U2709EN-BB
ELECTRONIC DISTRIBUTION