




The Multi-Mode AED for Professional Use


Because your patient's safety
is our number one priority, too.



Lifeline PRO



Manual



ECG



AED




defibtech

Life-Saving Technology Within Reach

Fast and simple



Manual Mode - Displays ECG and allows defibrillation energy selection (25-200J) with charge, shock, and disarm control.



ECG Mode - Displays ECG waveform using an optional 3-wire (LEAD II) ECG adaptor cable.



AED Mode - Provides fully-animated step-by-step video instructions for performing CPR, rescue breathing, and semi-automatic external defibrillation.



The Lifeline PRO AED is cleared by the FDA

